

For Research Use Only

SOHLH2 Polyclonal antibody

Catalog Number: 17057-1-AP



Basic Information

Catalog Number:

17057-1-AP

Size:

150ul, Concentration: 140 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10768

GenBank Accession Number:

BC025383

GeneID (NCBI):

54937

UNIPROT ID:

Q9NX45

Full Name:

spermatogenesis and oogenesis specific basic helix-loop-helix 2

Calculated MW:

225 aa, 25 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

SOHLH2, also named as TEB1, is a germ cell-specific bHLH transcription factor which may be involved in spermatogenesis and oogenesis. SOHLH2 has three isoforms with MW 25; 47 and 56 kDa.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

