

For Research Use Only

KATNAL2 Polyclonal antibody

Catalog Number: 17074-1-AP



Basic Information

Catalog Number: 17074-1-AP	GenBank Accession Number: BC034999	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 750 µg/ml by Nanodrop;	GeneID (NCBI): 83473	Recommended Dilutions: IF 1:200-1:800
Source: Rabbit	Full Name: katanin p60 subunit A-like 2	
Isotype: IgG	Calculated MW: 466 aa, 53 kDa	
Immunogen Catalog Number: AG10950	Observed MW: 40-80 kDa	

Applications

Tested Applications: IF, ELISA	Positive Controls: IF : NIH/3T3 cells,
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed NIH/3T3 cells using KATNAL2 antibody (17074-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).