

For Research Use Only

GMCL1L Polyclonal antibody

Catalog Number: 17129-1-AP



Basic Information

Catalog Number: 17129-1-AP	GenBank Accession Number: BC024185	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 107 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 64396	
Source: Rabbit	Full Name: germ cell-less homolog 1 (Drosophila)-like	
Isotype: IgG	Calculated MW: 526 aa, 60 kDa	
Immunogen Catalog Number: AG10873		

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

GMCL1L, also named as GCL and GMCL1P1, is possible function in spermatogenesis. It enhances the degradation of MDM2 and increases the amount of p53 probably by modulating the nucleocytoplasmic transport.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data