

For Research Use Only

ACTRT1 Polyclonal antibody, PBS Only

Catalog Number:17130-1-PBS



Basic Information

Catalog Number:

17130-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10874

GenBank Accession Number:

BC014597

GeneID (NCBI):

139741

UNIPROT ID:

Q8TDG2

Full Name:

actin-related protein T1

Calculated MW:

376 aa, 42 kDa

Observed MW:

40-45 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-ACTRT1 (IP:17130-1-AP, 4ug; Detection:17130-1-AP 1:800) with mouse testis tissue lysate 1800 ug. This data was developed using the same antibody clone with 17130-1-PBS in a different storage buffer formulation.