

For Research Use Only

# GALNTL5 Polyclonal antibody, PBS Only

Catalog Number: 17131-1-PBS



## Basic Information

Catalog Number:

17131-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10876

GenBank Accession Number:

BC022021

GeneID (NCBI):

168391

UNIPROT ID:

Q7Z4T8

Full Name:

UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylgalactosaminyltransferase-Like  
5

Calculated MW:

443 aa, 51 kDa

Observed MW:

51 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 17131-1-AP (GALNTL5 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17131-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human lung using 17131-1-AP (GALNTL5 antibody) at dilution of 1:100 (under 40x lens). This data was developed using the same antibody clone with 17131-1-PBS in a different storage buffer formulation.