

For Research Use Only

PLEKHG4 Polyclonal antibody

Catalog Number: 17156-1-AP **1 Publications**



Basic Information

Catalog Number: 17156-1-AP	GenBank Accession Number: BC082974	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 233 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 25894	
Source: Rabbit	Full Name: pleckstrin homology domain containing, family G (with RhoGef domain) member 4	
Isotype: IgG	Calculated MW: 1151 aa, 126 kDa	
Immunogen Catalog Number: AG10328		

Applications

Tested Applications:
ELISA

Cited Applications:
IF, IHC, WB

Species Specificity:
human, mouse

Cited Species:
human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Qingling Yuan	37336836	Apoptosis	WB,IHC,IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data