

For Research Use Only

# ARL17 Polyclonal antibody

Catalog Number: 17231-1-AP



## Basic Information

<b>Catalog Number:</b> 17231-1-AP	<b>GenBank Accession Number:</b> BC030570	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul	<b>GeneID (NCBI):</b> 641522	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q8IVW1	
<b>Isotype:</b> IgG	<b>Full Name:</b> ADP-ribosylation factor-like 17	
<b>Immunogen Catalog Number:</b> AG11072	<b>Calculated MW:</b> 160 aa, 17 kDa	

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data