

For Research Use Only

PUS7L Polyclonal antibody, PBS Only

Catalog Number:17273-1-PBS



Basic Information

Catalog Number:

17273-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11030

GenBank Accession Number:

BC033621

GeneID (NCBI):

83448

UNIPROT ID:

Q9H0K6

Full Name:

pseudouridylate synthase 7 homolog
(S. cerevisiae)-like

Calculated MW:

701 aa, 81 kDa

Observed MW:

81, 45 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17273-1-AP (PUS7L antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17273-1-PBS in a different storage buffer formulation.