

For Research Use Only

# LYSMD2 Polyclonal antibody

Catalog Number: 17292-1-AP



## Basic Information

Catalog Number:

17292-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11169

GenBank Accession Number:

BC022075

GeneID (NCBI):

256586

UNIPROT ID:

Q8IV50

Full Name:

LysM, putative peptidoglycan-binding, domain containing 2

Calculated MW:

124 aa, 14 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

