

For Research Use Only

ACTL8 Polyclonal antibody

Catalog Number: 17312-1-AP

1 Publications



Basic Information

Catalog Number:

17312-1-AP

Size:

150ul, Concentration: 500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11281

GenBank Accession Number:

BC028909

GeneID (NCBI):

81569

UNIPROT ID:

Q9H568

Full Name:

actin-like 8

Calculated MW:

366 aa, 41 kDa

Observed MW:

41 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : A549 cells, Jurkat cells, HeLa cells, HepG2 cells, MOLT-4 cells

Notable Publications

Author	Pubmed ID	Journal	Application
Xiaofei Zhang	28190767	Mol Cell	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17312-1-AP (ACTL8 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 17312-1-AP (ACTL8 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.