

For Research Use Only

ACTR8 Polyclonal antibody

Catalog Number: 17334-1-AP

1 Publications



Basic Information

Catalog Number:

17334-1-AP

Size:

150ul, Concentration: 270 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11367

GenBank Accession Number:

BC032744

GeneID (NCBI):

93973

UNIPROT ID:

Q9H981

Full Name:

ARP8 actin-related protein 8 homolog (yeast)

Calculated MW:

329 aa, 37 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

mouse

Positive Controls:

WB : A549 cells, Jurkat cells

Notable Publications

Author	Pubmed ID	Journal	Application
Jun Ho Lee	39243762	Cell	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17334-1-AP (ACTR8 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.