

For Research Use Only

TADA2B Polyclonal antibody, PBS Only

Catalog Number: 17367-1-PBS

Featured Product



Basic Information

Catalog Number:

17367-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10783

GenBank Accession Number:

BC013164

GeneID (NCBI):

93624

UNIPROT ID:

Q86TJ2

Full Name:

transcriptional adaptor 2 (ADA2 homolog, yeast)-beta

Calculated MW:

420 aa, 48 kDa

Observed MW:

55-62 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17367-1-AP (TADA2B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17367-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovary cancer tissue slide using 17367-1-AP (TADA2B antibody) at dilution of 1:400 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 17367-1-PBS in a different storage buffer formulation.