

For Research Use Only

GPBP1L1 Polyclonal antibody

Catalog Number: 17375-1-AP



Basic Information

Catalog Number: 17375-1-AP	GenBank Accession Number: BC022218	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 900 µg/ml by Nanodrop and 400 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 60313	Recommended Dilutions: WB 1:1000-1:5000
Source: Rabbit	Full Name: GC-rich promoter binding protein 1-like 1	
Isotype: IgG	Calculated MW: 474 aa, 52 kDa	
Immunogen Catalog Number: AG11180	Observed MW: 58 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HepG2 cells, Jurkat cells
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17375-1-AP (GPBP1L1 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.