

For Research Use Only

PSD2 Polyclonal antibody, PBS Only

Catalog Number: 17404-1-PBS

Featured Product



Basic Information

Catalog Number:

17404-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11386

GenBank Accession Number:

BC038233

GeneID (NCBI):

84249

UNIPROT ID:

Q9BQI7

Full Name:

pleckstrin and Sec7 domain containing 2

Calculated MW:

771 aa, 85 kDa

Observed MW:

85-95 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 17404-1-AP (PSD2 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17404-1-PBS in a different storage buffer formulation.