

For Research Use Only

TGIF2LX Polyclonal antibody

Catalog Number: 17405-1-AP



Basic Information

Catalog Number:

17405-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop and 300 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11424

GenBank Accession Number:

BC029920

GeneID (NCBI):

90316

UNIPROT ID:

Q8IUE1

Full Name:

TGFB-induced factor homeobox 2-like, X-linked

Calculated MW:

241 aa, 26 kDa

Observed MW:

34 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human testis tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 17405-1-AP (TGIF2LX antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.