

For Research Use Only

TGIF2LX Polyclonal antibody

Catalog Number: 17405-1-AP



Basic Information

Catalog Number: 17405-1-AP	GenBank Accession Number: BC029920	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 400 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 90316	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: TGFB-induced factor homeobox 2-like, X-linked	
Isotype: IgG	Calculated MW: 241 aa, 26 kDa	
Immunogen Catalog Number: AG11424	Observed MW: 34 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : human testis tissue,
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 17405-1-AP (TGIF2LX antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.