

For Research Use Only

CMBL Polyclonal antibody

Catalog Number: 17409-1-AP

1 Publications



Basic Information

Catalog Number:

17409-1-AP

Size:

150ul, Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11442

GenBank Accession Number:

BC001573

GeneID (NCBI):

134147

UNIPROT ID:

Q96DG6

Full Name:

carboxymethylenebutenolidase homolog (Pseudomonas)

Calculated MW:

245 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB: HeLa cells, mouse liver tissue, mouse skeletal muscle tissue, mouse kidney tissue, rat liver tissue

Notable Publications

Author	Pubmed ID	Journal	Application
Yingdan Huang	37967006	Cell Rep	WB, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17409-1-AP (CMBL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.