

For Research Use Only

BDNF Polyclonal antibody

Catalog Number: 17465-1-AP **7 Publications**



Basic Information

Catalog Number: 17465-1-AP	GenBank Accession Number: BC029795	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 307 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 627	
Source: Rabbit	UNIPROT ID: P23560	
Isotype: IgG	Full Name: brain-derived neurotrophic factor	
Immunogen Catalog Number: AG11329	Calculated MW: 247 aa, 28 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat

Cited Species:
rat, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Zhiyuan Wang	36438503	Theranostics	WB
Qiaozhen Qin	35883153	Stem Cell Res Ther	WB
Sen Yin	32802036	Neural Plast	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data