

For Research Use Only

BDNF Polyclonal antibody

Catalog Number: 17465-1-AP **13 Publications**



Basic Information

Catalog Number:

17465-1-AP

Size:

150ul, Concentration: 307 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11329

GenBank Accession Number:

BC029795

GeneID (NCBI):

627

UNIPROT ID:

P23560

Full Name:

brain-derived neurotrophic factor

Calculated MW:

247 aa, 28 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse, rat

Notable Publications

Author	Pubmed ID	Journal	Application
Zhiyuan Wang	36438503	Theranostics	WB
Qiaozhen Qin	35883153	Stem Cell Res Ther	WB
Sen Yin	32802036	Neural Plast	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

