

For Research Use Only

# ITGBL1 Polyclonal antibody

Catalog Number: 17484-1-AP

Featured Product

2 Publications



## Basic Information

Catalog Number:

17484-1-AP

Size:

150ul , Concentration: 193 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11521

GenBank Accession Number:

BC036788

GeneID (NCBI):

9358

UNIPROT ID:

O95965

Full Name:

integrin, beta-like 1 (with EGF-like repeat domains)

Calculated MW:

494 aa, 54 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

## Notable Publications

Author	Pubmed ID	Journal	Application
Zhe Wang	35596183	Biomed Eng Online	WB,IHC
Li Sun	27261588	Biomed Pharmacother	WB,IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

