

For Research Use Only

GNG11 Polyclonal antibody, PBS Only

Catalog Number:17503-1-PBS



Basic Information

Catalog Number: 17503-1-PBS	GenBank Accession Number: BC009709	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2791	
Source: Rabbit	UNIPROT ID: P61952	
Isotype: IgG	Full Name: guanine nucleotide binding protein (G protein), gamma 11	
Immunogen Catalog Number: AG11698	Calculated MW: 73 aa, 9 kDa	
	Observed MW: 9 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-11 (GBG11) is a subunit of Guanine nucleotide-binding proteins (G proteins). G proteins are involved as a modulator or transducer in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17503-1-AP (GNG11 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17503-1-PBS in a different storage buffer formulation.