

For Research Use Only

THNSL2 Polyclonal antibody

Catalog Number: 17573-1-AP



Basic Information

Catalog Number: 17573-1-AP	GenBank Accession Number: BC064423	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 800 µg/ml by Nanodrop;	GeneID (NCBI): 55258	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: threonine synthase-like 2 (S. cerevisiae)	
Isotype: IgG	Calculated MW: 405 aa, 45 kDa	
Immunogen Catalog Number: AG11788	Observed MW: 55-60 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HEK-293 cells, HepG2 cells
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17573-1-AP (THNSL2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.