

For Research Use Only

TSPYL4 Polyclonal antibody

Catalog Number: 17576-1-AP



Basic Information

Catalog Number: 17576-1-AP	GenBank Accession Number: BC009116	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 650 µg/ml by Nanodrop;	GeneID (NCBI): 23270	Recommended Dilutions: WB 1:2000-1:10000
Source: Rabbit	Full Name: TSPY-like 4	
Isotype: IgG	Calculated MW: 414 aa, 45 kDa	
Immunogen Catalog Number: AG10730	Observed MW: 45 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
mouse, rat

Positive Controls:

WB : mouse cerebellum tissue, rat cerebellum tissue

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17576-1-AP (TSPYL4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.