

For Research Use Only

GRHL1 Polyclonal antibody

Catalog Number: 17644-1-AP



Basic Information

Catalog Number:

17644-1-AP

Size:

150ul , Concentration: 200 ug/ml by Nanodrop and 167 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11832

GenBank Accession Number:

BC067519

GeneID (NCBI):

29841

UNIPROT ID:

Q9NZI5

Full Name:

grainyhead-like 1 (Drosophila)

Calculated MW:

618 aa, 70 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : NCCIT cell,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



NCCIT cell were subjected to SDS PAGE followed by western blot with 17644-1-AP (GRHL1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.