

For Research Use Only

PLA2G3 Polyclonal antibody, PBS Only

Catalog Number: 17646-1-PBS



Basic Information

Catalog Number: 17646-1-PBS	GenBank Accession Number: BC025316	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 50487	
Source: Rabbit	UNIPROT ID: Q9NZ20	
Isotype: IgG	Full Name: phospholipase A2, group III	
Immunogen Catalog Number: AG11845	Calculated MW: 509 aa, 57 kDa	

Applications

Tested Applications:
Indirect ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data