

For Research Use Only

# HBM Polyclonal antibody

Catalog Number: 17696-1-AP



---

## Basic Information

Catalog Number: 17696-1-AP	GenBank Accession Number: BC035682	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 3042	
Source: Rabbit	Full Name: hemoglobin, mu	
Isotype: IgG	Calculated MW: 141 aa, 16 kDa	
Immunogen Catalog Number: AG11948		

---

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human

---

## Background Information

---

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

---

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data