

For Research Use Only

MAN2B2 Polyclonal antibody, PBS Only

Catalog Number: 17697-1-PBS

Featured Product



Basic Information

Catalog Number:

17697-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11952

GenBank Accession Number:

BC033307

GeneID (NCBI):

23324

UNIPROT ID:

Q9Y2E5

Full Name:

mannosidase, alpha, class 2B, member 2

Calculated MW:

1009 aa, 114 kDa

Observed MW:

135 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 17697-1-AP (MAN2B2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17697-1-PBS in a different storage buffer formulation.



WB result of MAN2B2 antibody (17697-1-AP; 1:1200; incubated at room temperature for 1.5 hours) with sh-Control and sh-MAN2B2 transfected NIH/3T3 cells. This data was developed using the same antibody clone with 17697-1-PBS in a different storage buffer formulation.