

For Research Use Only

# VASN Polyclonal antibody

Catalog Number: 17714-1-AP



## Basic Information

Catalog Number: 17714-1-AP	GenBank Accession Number: BC068575	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 114990	
Source: Rabbit	Full Name: vasorin	
Isotype: IgG	Calculated MW: 673 aa, 72 kDa	
Immunogen Catalog Number: AG12040		

## Applications

Tested Applications:  
ELISA  
Species Specificity:  
human, mouse, rat

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data