

For Research Use Only

METTL4 Polyclonal antibody

Catalog Number: 17718-1-AP

2 Publications



Basic Information

Catalog Number: 17718-1-AP	GenBank Accession Number: BC111020	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 400 µg/ml by Nanodrop and 253 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 64863	
Source: Rabbit	Full Name: methyltransferase like 4	
Isotype: IgG	Calculated MW: 472 aa, 54 kDa	
Immunogen Catalog Number: AG12093		

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat

Cited Species:
human, mouse

Background Information

Methyltransferase like protein 4 (METTL4) is a internal m6Am methyltransferase, which mediates N6-methylation of Am30 on U2 snRNA in the context of an AAG motif in vivo and in vitro. METTL4 belongs to the MT-A70-like family, and like its paralogs METTL3/14, METTL4 is conserved from yeast to human.

Notable Publications

Author	Pubmed ID	Journal	Application
Xixi Tao	35781818	EMBO J	WB
Zhongmin Yang	33948366	Am J Cancer Res	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

