

For Research Use Only

ORC6 Polyclonal antibody, PBS Only

Catalog Number: 17784-1-PBS

Featured Product



Basic Information

Catalog Number:

17784-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11769

GenBank Accession Number:

BC063565

GeneID (NCBI):

23594

UNIPROT ID:

Q9Y5N6

Full Name:

origin recognition complex, subunit 6 like (yeast)

Calculated MW:

252 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17784-1-AP (ORC6 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17784-1-PBS in a different storage buffer formulation.