

For Research Use Only

LCN15 Polyclonal antibody

Catalog Number: 17793-1-AP



Basic Information

Catalog Number: 17793-1-AP	GenBank Accession Number: BC035124	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 350 ug/ml by Nanodrop;	GeneID (NCBI): 389812	Recommended Dilutions: IHC 1:50-1:500
Source: Rabbit	UNIPROT ID: Q6UWWO	
Isotype: IgG	Full Name: lipocalin 15	
Immunogen Catalog Number: AG11908	Calculated MW: 143 aa, 16 kDa	

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human stomach tissue,
Species Specificity: human	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

LCN15(lipocalin 15), a human-specific lipocalin family protein, as an abundant component of the olfactory mucus.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 17793-1-AP (LCN15 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).