

For Research Use Only

XAF1 Polyclonal antibody

Catalog Number: 17837-1-AP



Basic Information

Catalog Number: 17837-1-AP	GenBank Accession Number: BC073156	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 54739	
Source: Rabbit	Full Name: XIAP associated factor 1	
Isotype: IgG	Calculated MW: 301 aa, 35 kDa	
Immunogen Catalog Number: AG12110		

Applications

Tested Applications:
IF, ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of HeLa cells using 17837-1-AP (XAF1 antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).