

For Research Use Only

NEUROG3 Polyclonal antibody

Catalog Number: 17902-1-AP

1 Publications



Basic Information

Catalog Number:

17902-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12272

GenBank Accession Number:

BC074938

GeneID (NCBI):

50674

UNIPROT ID:

Q9Y4Z2

Full Name:

neurogenin 3

Calculated MW:

23 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Priyadharsini Raju	25661254	Reproduction	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

