

For Research Use Only

PCDH17 Polyclonal antibody

Catalog Number: 17915-1-AP



Basic Information

Catalog Number:

17915-1-AP

Size:

150ul , Concentration: 900 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12322

GenBank Accession Number:

BC028165

GeneID (NCBI):

27253

UNIPROT ID:

O14917

Full Name:

protocadherin 17

Calculated MW:

1159 aa, 126 kDa

Observed MW:

160-170 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17915-1-AP (PCDH17 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.