

For Research Use Only

SOX4 Polyclonal antibody

Catalog Number: 17919-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:
17919-1-AP

Size:
150ul, Concentration: 367 µg/ml by Bradford method using BSA as the standard;

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
AG12333

GenBank Accession Number:
BC072668

GeneID (NCBI):
6659

Full Name:
SRY (sex determining region Y)-box 4

Calculated MW:
474 aa, 47 kDa

Purification Method:
Antigen affinity purification

Applications

Tested Applications:
ELISA

Cited Applications:
IHC, WB

Species Specificity:
human, mouse, rat

Cited Species:
human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhen-Zhen Chen	27553854	Nat Commun	WB, IHC
Rongfei Han	27446439	Oncol Lett	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data