

For Research Use Only

SPRY3 Polyclonal antibody, PBS Only

Catalog Number:17932-1-PBS



Basic Information

Catalog Number:

17932-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12364

GenBank Accession Number:

BC074973

GeneID (NCBI):

10251

UNIPROT ID:

O43610

Full Name:

sprouty homolog 3 (Drosophila)

Calculated MW:

288 aa, 31 kDa

Observed MW:

31-33 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



RAW264.7 cells were subjected to SDS PAGE followed by western blot with 17932-1-AP (SPRY3 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17932-1-PBS in a different storage buffer formulation.