

For Research Use Only

# GABRP Polyclonal antibody

Catalog Number: 17957-1-AP



## Basic Information

Catalog Number:	17957-1-AP	GenBank Accession Number:	BC074865	Purification Method:	Antigen affinity purification
Size:	150ul	GeneID (NCBI):	2568		
Source:	Rabbit	Full Name:	gamma-aminobutyric acid (GABA) A receptor, pi		
Isotype:	IgG	Calculated MW:	440 aa, 51 kDa		
Immunogen Catalog Number:	AG12453				

## Applications

Tested Applications:  
ELISA  
Species Specificity:  
human, mouse, rat

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1%BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data