

For Research Use Only

ADH1A Polyclonal antibody

Catalog Number: 17998-1-AP

1 Publications



Basic Information

Catalog Number:

17998-1-AP

Size:

150ul, Concentration: 350 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12467

GenBank Accession Number:

BC074738

GeneID (NCBI):

124

UNIPROT ID:

P07327

Full Name:

alcohol dehydrogenase 1A (class I), alpha polypeptide

Calculated MW:

375 aa, 40 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:12000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse, rat

Positive Controls:

WB: mouse liver tissue, mouse lung tissue, rat liver tissue, rat lung tissue

Notable Publications

Author	Pubmed ID	Journal	Application
Zihao Jiang	39914624	Toxicol Appl Pharmacol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17998-1-AP (ADH1A antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.