

For Research Use Only

APOA5 Polyclonal antibody, PBS Only

Catalog Number:18019-1-PBS



Basic Information

Catalog Number:

18019-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12568

GenBank Accession Number:

BC101787

GeneID (NCBI):

116519

UNIPROT ID:

Q6Q788

Full Name:

apolipoprotein A-V

Calculated MW:

366 aa, 41 kDa

Observed MW:

42-45 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cell lysates were subjected to SDS PAGE followed by western blot with 18019-1-AP (APOA5 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18019-1-PBS in a different storage buffer formulation.