

For Research Use Only

# RPS6KC1 Polyclonal antibody

Catalog Number: 18020-1-AP



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
18020-1-AP	BC104769	Antigen affinity purification
Size:	GeneID (NCBI):	
150ul	26750	
Source:	Full Name:	
Rabbit	ribosomal protein S6 kinase, 52kDa, polypeptide 1	
Isotype:	Calculated MW:	
IgG	1066 aa, 119 kDa	
Immunogen Catalog Number:		
AG12569		

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human, mouse, rat

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data