

For Research Use Only

GREM1/Gremlin 1 Polyclonal antibody, PBS Only

Catalog Number: 18024-1-PBS



Basic Information

Catalog Number:

18024-1-PBS

Size:

100ug, Concentration: 1mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12588

GenBank Accession Number:

BC093778

GeneID (NCBI):

26585

UNIPROT ID:

O60565

Full Name:

gremlin 1, cysteine knot superfamily,
homolog (Xenopus laevis)

Calculated MW:

184 aa, 21 kDa

Observed MW:

21 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 18024-1-AP (GREM1/Gremlin 1 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18024-1-PBS in a different storage buffer formulation.