

For Research Use Only

KCNK9 Polyclonal antibody, PBS Only

Catalog Number:18033-1-PBS



Basic Information

Catalog Number:

18033-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12603

GenBank Accession Number:

BC075079

GeneID (NCBI):

51305

UNIPROT ID:

Q9NPC2

Full Name:

potassium channel, subfamily K, member 9

Calculated MW:

374 aa, 42 kDa

Observed MW:

45-50 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 18033-1-AP (KCNK9 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18033-1-PBS in a different storage buffer formulation.