For Research Use Only

PCDHB3 Polyclonal antibody

Catalog Number: 18076-1-AP



Basic Information

Catalog Number:

GenBank Accession Number:

Purification Method: Antigen affinity purification

18076-1-AP

Size:

BC104817 GeneID (NCBI): 56132

150ul Source:

UNIPROT ID: Q9Y5E6

Rabbit Isotype:

Full Name: protocadherin beta 3

Immunogen Catalog Number:

Calculated MW:

796 aa, 87 kDa

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Storage

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

Selected Validation Data