

For Research Use Only

HERV-FRD Polyclonal antibody, PBS Only

Catalog Number: 18101-1-PBS



Basic Information

Catalog Number:

18101-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12380

GenBank Accession Number:

BC068585

GeneID (NCBI):

405754

UNIPROT ID:

P60508

Full Name:

HERV-FRD provirus ancestral Env polyprotein

Calculated MW:

538 aa, 60 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 18101-1-AP (HERV-FRD antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18101-1-PBS in a different storage buffer formulation.