

For Research Use Only

BCL2L10 Polyclonal antibody, PBS Only

Catalog Number: 18114-1-PBS



Basic Information

Catalog Number:

18114-1-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12596

GenBank Accession Number:

BC093826

GeneID (NCBI):

10017

UNIPROT ID:

Q9HD36

Full Name:

BCL2-like 10 (apoptosis facilitator)

Calculated MW:

204 aa, 23 kDa

Observed MW:

22 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data