

For Research Use Only

# CNO Polyclonal antibody

Catalog Number: 18217-1-AP



---

## Basic Information

Catalog Number:  
18217-1-AP

Size:  
150ul

Source:  
Rabbit

Isotype:  
IgG

Immunogen Catalog Number:  
AG12975

GenBank Accession Number:  
BC067815

GeneID (NCBI):  
55330

Full Name:  
cappuccino homolog (mouse)

Calculated MW:  
217 aa, 23 kDa

Purification Method:  
Antigen affinity purification

---

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human, mouse, rat

---

## Background Information

---

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

---

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data