

For Research Use Only

GADL1 Polyclonal antibody, PBS Only

Catalog Number: 18240-1-PBS



Basic Information

Catalog Number:

18240-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13005

GenBank Accession Number:

BC093701

GeneID (NCBI):

339896

UNIPROT ID:

Q6ZQY3

Full Name:

glutamate decarboxylase-like 1

Calculated MW:

39 kDa, 59 kDa, 47 kDa

Observed MW:

47 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with 18240-1-AP (GADL1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18240-1-PBS in a different storage buffer formulation.



IP result of anti-GADL1 (IP:18240-1-AP, 4ug; Detection:18240-1-AP 1:300) with mouse brain tissue lysate 4000ug. This data was developed using the same antibody clone with 18240-1-PBS in a different storage buffer formulation.