

For Research Use Only

# NFATC1 Polyclonal antibody

Catalog Number: 18337-1-AP **3 Publications**



## Basic Information

Catalog Number:

18337-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13014

GenBank Accession Number:

BC104753

GeneID (NCBI):

4772

UNIPROT ID:

O95644

Full Name:

nuclear factor of activated T-cells,  
cytoplasmic, calcineurin-dependent 1

Calculated MW:

716aa,78 kDa; 943aa,101 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

## Notable Publications

Author	Pubmed ID	Journal	Application
Lin Luo	37783002	Biomater Adv	WB
Ningning Jiang	37473680	Biomed Pharmacother	WB
Wengang Wang	37408034	Clin Epigenetics	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data