

For Research Use Only

NME4 Polyclonal antibody

Catalog Number: 18348-1-AP



Basic Information

Catalog Number: 18348-1-AP	GenBank Accession Number: BC004880	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 4833	
Source: Rabbit	UNIPROT ID: O00746	
Isotype: IgG	Full Name: non-metastatic cells 4, protein expressed in	
Immunogen Catalog Number: AG13162	Calculated MW: 21 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Notable Publications

Author	Pubmed ID	Journal	Application
Lu Huang	25463732	Environ Pollut	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

