

For Research Use Only

CYBA Polyclonal antibody

Catalog Number: 18440-1-AP



Basic Information

Catalog Number:

18440-1-AP

Size:

150ul, Concentration: 1500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13071

GenBank Accession Number:

BC006465

GeneID (NCBI):

1535

UNIPROT ID:

P13498

Full Name:

cytochrome b-245, alpha polypeptide

Calculated MW:

21 kDa

Observed MW:

22 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : 37°C incubated A549 cells, 37°C incubated HepG2 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 18440-1-AP (CYBA antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.